



All of the biometric images and graphs produced by cardiological diagnostics (angiographs, polygraphs, echocardiographs, EKG, Holter, etc.) are acquired and saved in the DicomWare database, viewable using Medweb, and reported using eVisit: a unique solution, covering the whole cardiology field.

Elemento	Data	Ora	Descrizione	Evento
Evento	11/03/2010	09:39	ingresso in sala	Accesso
Accesso	11/03/2010	09:40	percutaneo	Farmaco
Farmaco	11/03/2010	09:40	Eparina 12 ml	Procedura
Farmaco	11/03/2010	09:40	Atenololo 2 ml	Materiale
Materiale	11/03/2010	09:41	xxx	Lesione
Lesione	11/03/2010	09:42	ristenosi su BMS su Coronaria destra al tratto medio [CDx med]	Complicanze
Lesione	11/03/2010	09:43	ristenosi su DES su arteria mammaria (su ...)	Par. vitali
REFERTO	11/03/2010	09:43		IMA in atto

eVisit can import all of the measurements taken with **echocardiographs, angiographs, EKG's** directly on-screen in **Structured DICOM Report** format.

eVisit also provides an environment for assessing the main needs in echocardiography, hemodynamics and electrocardiography, and can automatically fill out the report forms with the results shown on the image and graph viewer.

Lesione **ristenosi su BMS su Coronaria destra al tratto medio [CDx med]**

Lesione **Device**

Data e ora: 11/03/2010 09:42 Numero lesione: 1

Ramo: Vaso nativo

Segmento: Coronaria destra al tratto medio [CDx med]

Tipo lesione: ristenosi su BMS Classe ristenosi: III

Ostiale: si

Lesione biforcata:

Colpevole di SCA: no

Ulcerazione: si

Trombo: no

Calificazioni: essenti

Tortuosita': marcata

Angolazione: < 45°


Occlusione: recente

Ultimo vaso residuo: si

No Reflow: si

**DIMENSIONI**

BASALE	FINALE
Lunghezza (mm)	Lunghezza (mm)
Diametro (mm)	Diametro (mm)
Percentuale	Percentuale
Flusso TIMI	Flusso TIMI



Displaying EKG and CR chest images for pre-operative assessment